

# Bleed test of septa by UV-Visible


## Experimental goal

Bleed test of septa in the ACN/H<sub>2</sub>O solvent by UV-Visible detection.


EXPERIMENTAL	
instrument	EVOLUTION 220,UV-Visble Spectrophotometer
Wavelength	190 nm-500 nm
Integration time	0.05 sec
Solvent	ACN:H <sub>2</sub> O=1:1

## Experimental result

Blank:ACN:H<sub>2</sub>O=1:1




Finetech



<b>nm(Finetech) 1</b>	<b>Abs</b>
203.86	1.650
296.92	0.388
<b>nm(Finetech) 2</b>	<b>Abs</b>
203.92	1.648
296.92	0.389
<b>nm(Finetech) 3</b>	<b>Abs</b>
203.96	1.650
296.92	0.396
<b>nm(Finetech)</b>	<b>RSD%</b>
203.9	0.070
296.9	1.115

## Brand A



nm( BrandA )	Abs
205.61	1.586
296.97	0.377

## Summary

After 24 hours, the signal intensity of salicylic acids was not interference by extract from both septa. The RSD% of signal intensity for both septa were less than 1%. Use of this lower bleed materials reduces the potential for sample error caused by septa bleed and improve overall reliability in the ACN :H<sub>2</sub>O(50:50) solvent condition.

nm(Finotech) 1	Abs	nm(BrandA)	Abs
203.86	1.650	205.61	1.586
296.92	0.388	296.97	0.377
nm(Finotech) 2	Abs		
203.92	1.648		
296.92	0.389		
nm(Finotech) 3	Abs		
203.96	1.650		
296.92	0.396		
nm(Finotech)	RSD%		
203.9	0.070		
296.9	1.115		