

fine ech FINETECH RESEARCH & INNOVATION CORP.

Cell Culture Series



ABOUT US

Constantly researching and innovating...

Finetech Research and Innovation Corporation was established in 1999 and has been engaged in the gas and fluid filtration used in laboratories and the medical fields for many years. We concentrate on the researching, developing, and manufacturing of lab filters such as syringe filters, 50mm air/vent filters, pipette filters, gas analyzer filters, and customized filters.

In 2004, we were the first manufacturer of transducer protectors used in dialysis. In 2017, sales of transducer protectors had reached to over 30 million pieces/year. Our high quality and ISO 13485 certified transducer protectors has kept our customers satisfied for more than 14 years.

Manufactured in 2010, our syringe filter is designed to provide efficient filtration of aqueous and organic solutions. It can be used to remove particles from a sample prior to analysis by HPLC, for the filtration of gases, and for the removal of bacteria from a sample. Syringe filters are widely used in environmental, pharmaceutical, biotechnology, and agricultural testing laboratories.

In 2015, Finetech designed and manufactured caps, septa, and samples vials used in chromatography. We supply HPLC consumables to the global chromatography community in over 36 countries. We constantly innovate with releases of new lab consumables and medical devices every year. Some of the other products we manufacture are sterile syringe filter, membrane filters, centrifuge tube, glass microfiber filters, vacuum filtration system, and QuECh-ERS. To find out more, please visit our website "www.finetech-filters.com".

History

- 2021 Developed Racked Filters Tips, PCR Supplies
- **2020 Produced Disposable Medical Masks**
- 2019 Produced Bottle Top Filter Unit
- 2018 Researched and developed QuEChERS
- 2017 Produced 15mL & 50mL centrifuge tubes (sterile & non-sterile)
- 2016 New product launched sterile filtration cup
- 2015 Achieved ISO 9001. Received award of Conventional Industry
 Technology Development for Medical grade septa
- 2014 Started producing sterile syringe filters
- 2012 Produced the 1st 50mm air vent filters
- 2010 Produced the 1st syringe filters manufactured in a cleanroom
- 2008 Developed transducer protectors. 1st manufacturer in Asia
- 2005 Awarded the Medical Device award from Taiwan Government
- 1999 Finetech Research and Innovation Corporation established









Taiwan Factory



China Office




Vietnam 1st factory



New factory under construction in Vietnam



Product Timeline



Contents

A. Finetech® Sterile Syringe filter04
B. Finetech® Bottle Top Filter05
C. Finetech® Cell Culture Plate06
D. Finetech® Cell Culture Flask07
E. Finetech® Cell Culture Dish08
F. Finetech® Sterile Centrifuge Tube09
G. Finetech Pipette Filter Tips11
H. Finetech® Sterile Serological Pipette13



A. Finetech[®] Sterile Syringe filter



• Housing material: Medical grade polypropylene

Pore Size: 0.22um, 0.45µm
Outside diameter: 17mm
Membrane diameter: 13mm

• High flow rate and throughput

• Low hold-up volume (less than 25μL)

Low extractables

• Individual sterile blister pack



• Housing material: Medical grade polypropylene

Pore Size: 0.22um, 0.45µm
Outside diameter: 33mm
Membrane diameter: 25mm
High flow rate and throughput

• Low hold-up volume (less than 100µL)

Low extractables

• Individual sterile blister pack

Ordering Information

① Membrane	② Diameter	③ Sterile	④ Pore size	(5) Membrane characteristics
PVDF CA PTFE PES Nylon RC MCE GF	017=13mm	N (Non-Sterile)	022=0.22μm	I = hydrophilic
	033=33mm	S (EO/GAMMA)	045=0.45μm	O = hydrophobic

13mm: 100pcs/pk 33mm: 100pcs/pk

Example: PVDF $033 \le 022 I = 33$ mm, PVDF membrane, Sterile, 0.22um

04

B. Finetech® Bottle Top Filter



Features

- Manufactured from virgin Polystyrene.
- Tissue culture treated surface to aid in cell attachment.
- Walls have even thickness and no distortion.
- Wide design increases stability when stacked.
- Quick opening cap.
- Angled bottleneck, convenient for pipettes and cell scraping.
- Frosted write-on spot.
- Clear graduation scale.



	Membrane	Diameter	Pore Size	Description	Packaging
FC-PVDF500G022I-90	PVDF	90mm	0.22um	Sterile Bottle Top Vacuum Filter Funnel Capacity: 500ml Individually Packed	12pcs/CTN
FC-PVDF500G0451-90	PVDF	90mm	0.45um	Sterile Bottle Top Vacuum Filter Funnel Capacity: 500ml Individually Packed	12pcs/CTN
FC-PES500G022I-90	PES	90mm	0.22um	Sterile Bottle Top Vacuum Filter Funnel Capacity: 500ml Individually Packed	12pcs/CTN
FC-PES500G0451-90	PES	90mm	0.45um	Sterile Bottle Top Vacuum Filter Funnel Capacity: 500ml Individually Packed	12pcs/CTN

C. Finetech® Cell Culture Plate





Features

- Available in 6, 12, 24, 48, 96 well plates.
- Manufactured from virgin Polystyrene.
- Tissue culture treated surface to aid in cell attachment.
- Tight fitting cover prevents contamination and loss by evaporation.
- One-way cover reduces errors.
- Alphanumeric labels for easy organization.
- Stack-able to save space.
- E-Beam sterilization.
- DNase/RNase-free and non-pyrogenic.

Model No.	Number of Wells	Bottom	Cell Growth Area	Recommended Volume	Sterile	Packaging
CCWP6-S	6	Flat	9.5 cm²	1.9-2.9 mL	Yes	50pcs/carton
CCWP12-S	12	Flat	3.6 cm²	0.76-1.14 mL	Yes	50pcs/carton
CCWP24-S	24	Flat	1.9 cm²	0.38-0.57 mL	Yes	50pcs/carton
CCWP48-S	48	Flat	0.88 cm ²	0.19-0.285 mL	Yes	50pcs/carton
CCWP96-S	96	Flat	0.32 cm ²	0.1-0.2 mL	Yes	65pcs/carton

D. Finetech[®] Cell Culture Flask






Features

- Manufactured from virgin Polystyrene.
- Tissue culture treated surface to aid in cell attachment.
- Walls have even thickness and no distortion.
- Wide design increases stability when stacked.
- Quick opening cap.
- Angled bottleneck, convenient for pipettes and cell scraping.
- Frosted write-on spot.
- Clear graduation scale.
- Stack-able to save space.
- E-Beam sterilization.

Model No.	Volume	Recommende d Volume	Cell Growth Area	Cap	Sterile	Packaging
CCF25SC-S	50 mL	5-7.5 mL	25 cm²	Sealed	Yes	12pcs/pack 300pcs/case
CCF75SC-S	250 mL	15-22.5 mL	75 cm²	Sealed	Yes	5pcs/pack 90pcs/case
CCF175SC-S	750 mL	35-52.5 mL	175 cm²	Sealed	Yes	5pcs/pack 50pcs/case
CCF25FC-S	50 mL	5-7.5 mL	25 cm²	Vent	Yes	12pcs/pack 300pcs/case
CCF75FC-S	250 mL	15-22.5 mL	75 cm²	Vent	Yes	5pcs/pack 90pcs/case
CCF175FC-S	750 mL	35-52.5 mL	175 cm²	Vent	Yes	5pcs/pack 50pcs/case

E. Finetech[®] Cell Culture Dish





Features

- Tissue culture treated surface to aid in cell attachment.
- Flat and with even thickness. Allowing distortion-free observation.
- Easy to grip sides.
- E-Beam sterilization.
- DNase/RNase-free and non-pyrogenic.

Model No.	Diameter	Recommended Volume	Cell Growth Area	Sterile	Packaging
CCD35-S	35mm	1.8-2.7 mL	8.5 cm²	Yes	20pcs/pack
	33,,,,,,		0.0 0	. 00	500pcs/case
000/0.0	(0)	4.2.7.2 mal	22.9 cm ²	Vaa	20pcs/pack
CCD60-S	60mm	4.2-6.3 mL	22.9 CIII-	Yes	500pcs/case
000100 0	100mm	44.47.5	F7 /2	V	10pcs/pack
CCD100-S	(Grip Design)	11-16.5 mL	57.6 cm ²	Yes	300pcs/case

F. Finetech® Sterile Centrifuge Tube

15/50ml Centrifuge Tubes




Features


- These conical-bottom tubes are made of ultra-clear polypropylene.
- Cap Material: HDPE
- Sterility: EO or Gamma
- Leakproof
- Black printed graduations with large marking strip for easy labeling.
- Easy single-hand operation.
- Clear graduations provides easy volume reference.
- Available in non-sterile.

	15mL	50mL
Bottom Graduation Range	0.1 ~ 1 mL	-
Bottom Graduation Interval	0.1 mL	-
Graduation Range	1.5 ~ 14.5 mL	5 ~ 50 mL
Graduation Interval	0.5 mL	2.5 mL
Freezable	-80°C	-80°C
Maximum RCF	21,800 xg	14,000 xg
Maximum RPM	15,000 rpm	12,000 rpm

Model No.	Description	Packaging
ST-tube15	15 mL centrifuge tubes, Nonsterile	25pcs/pk, 500pcs/ctn
ST-tube50	50 mL centrifuge tubes, Nonsterile	25pcs/pk, 500pcs/ctn
ST-tube15-EO	15 mL centrifuge tubes, EO sterile	25pcs/pk, 500pcs/ctn
ST-tube50-EO	50 mL centrifuge tubes, EO sterile	25pcs/pk, 500pcs/ctn
ST-tube15-GA	15 mL centrifuge tubes, Gamma sterile	25pcs/pk, 500pcs/ctn
ST-tube50-GA	50 mL centrifuge tubes, Gamma sterile	25pcs/pk, 500pcs/ctn



G. Finetech[®] Pipette Filter Tips



Features

- Made of highly transparent PP material.
- Non-pyrogenic and with low adhesion.
- Compatible with Finetech, Eppendorf, Brand, Gilson, Thermo/Finn, and other brands pipette.
- High quality tip with smooth inner walls reduces sample residual.
- Racked sterile filter tips also available.
- Free of DNase, RNase, DNA, Endotoxin.

Racked pipette tips with filter







Model No.	Description	Packaging
Tip-S10F-R	10ul Sterile filter pipette tips, Clear, Racked	96/pk,4800/cs
Tip-S10FL-R	10ul Sterile filter pipette tips, Low Retention, Clear, Racked	96/pk,4800/cs
Tip-S10FE-R	10ul Sterile Extended Long pipette tips, Clear, Racked	96/pk,4800/cs
Tip-S10FEL-R	10ul Sterile Extended Long pipette tips, Low Retention, Clear, Racked	96/pk,4800/cs
Tip-S200F-R	200ul Sterile filter pipette tips, Clear, Racked	96/pk,4800/cs
Tip-S200FL-R	200ul Sterile filter pipette tips, Low Retention, Clear, Racked,	96/pk,4800/cs
Tip-S300F-R	300ul Sterile filter pipette tips, Clear, Racked	96/pk,4800/cs
Tip-S300FL-R	300ul Sterile filter pipette tips, Low Retention, Clear, Racked,	96/pk,4800/cs
Tip-S1000F-R	1000ul Sterile filter pipette tips, blue, Racked	96/pk,4800/cs
Tip-S1000FL-R	1000ul Sterile filter pipette tips, Low Retention, Clear, Racked	96/pk,4800/cs
Tip-S1250F-R	1250ul Sterile filter pipette tips, blue, Racked	96/pk,4800/cs
Tip-S1250FL-R	1250ul Sterile filter pipette tips, Low Retention, Clear, Racked	96/pk,4800/cs


H. Finetech® Sterile Serological Pipette



- Capacity:1mL-50mL
- Sterilization:E-Beam
- Packaging:Individually packed

Features

- made of imported general purpose polystyrene (GPPS) in clean room environment.
- Available in six capacity: 1.0ml, 2.0ml, 5.0ml, 10.0ml, 25.0ml, 50.0ml. Graduations are calibrated within ±2%.
- 1.0ml, 2.0ml, 5.0ml pipettes are shaped by streching method.
- 5.0ml, 10.0ml, 25.0ml, 50.0ml pipettes tips are welded to the tube by supersonic wave.
- Smooth and clean without draw lines.
- Different colors correspond to different specification, making it easy to distinguish.
- Maximum reduction of liquid adhesion to the tube inner wall, which improve the sampling accuracy.
- Excellent double-side graduation design ideally meet with the demand for liquid adding or reducing.
- Clear and simple graduation, as well as the negative one for extra volume, makes the liquid drawing and reading easier.
- All serological pipettes are assembled with high quality filters.
- Available sterile or non-sterile .
- Non-pyrogenic and DNase/RNase-free.



Model No.	Description	Packaging
	Capacity: 1mL	
FT9010	Size: Ø5×270mm	
	Individual packs and E-Beam sterilization	300pcs/box ,1800pcs/CTN
	Capacity: 2mL	300pcs/box , 1800pcs/CTN
FT9020	Size: Ø6×280mm	
	Individual packs and E-Beam sterilization	
	Capacity: 5mL	
FT9050	Size: Ø8×345mm	
	Individual packs and E-Beam sterilization	200mas /hav. 1200mas /CTN
	Capacity: 10mL	200pcs/box 1200pcs/CTN
FT9100	Size: Ø10×345mm	
	Individual packs and E-Beam sterilization	
	Capacity: 25mL	
FT9250	Size: Ø15×345mm	100pcs/box 600pcs/CTN
	Individual packs and E-Beam sterilization	
	Capacity: 50mL	
FT9500	Size: Ø18×345mm	75pcs/box 450pcs/CTN
	Individual packs and E-Beam sterilization	

FINETECH RESEARCH & INNOVATION CORP.



Company Profile

- Since 1999, Finetech Research and Innovation Corporation has been manufacturing and supplying laboratory consumables to various types of labs all over the world. With exports to over 40 countries, Finetech has become one of the leading provider for environmental, pharmaceutical, biological, and other testing laboratories
- Finetech began with the research & development and manufacturing of syringe filters and other filtration products. In keeping with our core principle of innovation, Finetech has expanded their product line to include lab equipment and other supplies to meet the needs of their customers.
- Finetech follows the principles of quality over quantity and good services for all.

The importance of quality can be seen in the use of quality imported materials and equipment, the use of a cleanroom, and the strict quality control standards.







Contact us

FINETECH RESEARCH & INNOVATION CORP



